

Ucp3 Rabbit Pab For Wb, Ihc-P, Elisa - P55916

Item Number: CM0019577

Introduction

Rabbit polyclonal antibody against human UCP3, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Suitable for studying mitochondrial metabolism and uncoupling.

[Learn More](#)

Parameter	Specification
Product Name	UCP3 Rabbit pAb
Remarks/Alias	SLC25A9; UCP3
Species	Human
GeneID (Human)	7352
GeneID	7352
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-119 of human UCP3 (NP_073714.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55916
Protein Weight	34kDa
Shipping	Ice bag