


Ufd1L Rabbit Mab - Q92890

Item Number: CM0003826



Introduction

Rabbit monoclonal antibody against UFD1L (ubiquitin recognition factor in ERAD) validated for Western blot and ELISA in human and rat samples. Ideal for studying ubiquitin-proteasome pathways.

[Learn More](#)

Parameter	Value
Product Name	UFD1L Rabbit mAb
Remarks/Alias	UFD1L
Species	Human
Gene ID (Human)	7353
Gene ID	7353
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q92890
Protein Weight	35kDa
Shipping	Ice bag