

Ugt2B7 Rabbit Pab - P16662

Item Number: CM0022216



Introduction

Rabbit polyclonal antibody against UGT2B7. Validated in WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets amino acids 25-351. Ideal for drug metabolism and hormone conjugation research.

[Learn More](#)

Parameter	Value
Product Name	UGT2B7 Rabbit pAb
Remarks/Alias	UGT2B9; UDPGTH2; UDPGT2B7; UDPGTh-2; UDPGT 2B7; UDPGT 2B9; UGT2B7
Species	Human
GeneID (Human)	7364
GeneID	7364
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-351 of human UGT2B7 (NP_001065.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16662
Protein Weight	61kDa
Shipping	Ice bag