

Ulbp2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bzm5

Item Number: CM0025684



Introduction

ULBP2 rabbit polyclonal antibody for Western blot and ELISA. Recognizes human ULBP2/NKG2D ligand 2, expressed in cancer and fetal tissues, mediating NK cell activation. Ideal for immuno-oncology studies.

[Learn More](#)

Parameter	Value
Product Name	ULBP2 Rabbit pAb
Remarks/Alias	N2DL2, RAET1H, RAET1L, NKG2DL2, ALCAN-alpha, ULBP2
Species	Human
Gene ID	80328
Immunogen	A synthetic peptide corresponding to amino acids 100-200 of human ULBP2 (NP_079493.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BZM5
Protein Weight	27kDa
Shipping	Ice bag