

## Unc5C Rabbit Polyclonal Antibody For Wb And Elisa - 095185

Item Number: CM0013449

### Introduction

High-quality UNC5C rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat UNC5C (103 kDa). Ideal for axon guidance and apoptosis studies.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | UNC5C Rabbit pAb                                                                                                   |
| Remarks/Alias    | UNC5H3; UNC5C                                                                                                      |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 8633                                                                                                               |
| GeneID           | 8633                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 300-360 of human UNC5C (NP_003719.3) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | O95185                                                                                                             |
| Protein Weight   | 103kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                            |