

Upf1/Rent1 Rabbit Pab - Q92900

Item Number: CM0015188



Introduction

Rabbit polyclonal antibody against human UPF1/RENT1, validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Suitable for nonsense-mediated mRNA decay (NMD) research.

[Learn More](#)

Parameter	Value
Product Name	UPF1/RENT1 Rabbit pAb
Remarks/Alias	UTF; HUPF1; NORF1; RENT1; smg-2; pNORF1; UPF1/RENT1
Species	Human
GeneID (Human)	5976
GeneID	5976
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 819-1118 of human UPF1/RENT1 (NP_002902.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92900
Protein Weight	124 kDa
Shipping	Ice bag