


Uqcrh Rabbit Mab For Wb, Ihc-P, Elisa - P07919

Item Number: CM0010274



Introduction

Recombinant rabbit monoclonal antibody targeting human UQCRH, a mitochondrial complex III subunit. Validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for mitochondrial disease and bioenergetics research.

[Learn More](#)

Parameter	Value
Product Name	UQCRH Rabbit mAb
Remarks/Alias	QCR6; UQCR8; MC3DN11; UQCRH
Species	Human
GeneID (Human)	7388
GeneID	7388
Immunogen	Synthetic peptide corresponding to amino acids 1-91 of human UQCRH (P07919)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07919
Protein Weight	11kDa
Shipping	Ice bag