


Uromodulin/Umod Rabbit Mab For Wb, If-P, Ihc-P, Elisa - P07911

Item Number: CM0006981



Introduction

Recombinant rabbit monoclonal antibody targeting human Uromodulin/UMOD (THP), validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with human and mouse. Ideal for renal biology, UTI research, and biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	Uromodulin/UMOD Rabbit mAb
Remarks/Alias	THP; FJHN; HNFJ; THGP; HNFJ1; MCKD2; ADTKD1; ADMCKD2; Uromodulin/UMOD
Species	Human
GeneID (Human)	7369
GeneID	7369
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 175-324 of human Tamm Horsfall Glycoprotein(NP_003352.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P07911
Protein Weight	70kDa/55kDa/115kDa
Shipping	Ice bag