

## Ush1C Rabbit Pab - Ush1C

Item Number: CM0024824

### Introduction

Polyclonal antibody targeting zebrafish USH1C, validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Suitable for research on Usher syndrome and sensory neuroscience.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | USH1C Rabbit pAb                                                                                                         |
| Remarks/Alias    | zgc:136806; USH1C                                                                                                        |
| Species          | Zebrafish                                                                                                                |
| Gene ID (Human)  | 10083                                                                                                                    |
| Gene ID          | 564412                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of zebrafish USH1C (NP_001035018.1). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, ELISA                                                                                                                |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| Protein Weight   | 45kDa/58kDa/60kDa/62kDa/101kDa                                                                                           |
| Shipping         | Ice bag                                                                                                                  |