

Usp10 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q14694

Item Number: CM0008843



Introduction

High-quality rabbit monoclonal antibody targeting human USP10 (Q14694). Validated for WB, IF/ICC, and ELISA applications. Cross-reacts with mouse. Ideal for studying deubiquitination, autophagy, and DNA damage response pathways.

[Learn More](#)

Parameter	Value
Product Name	USP10 Rabbit mAb
Remarks/Alias	UBPO; UBP10; MGC2621; USP10
Species	Human
Gene ID (Human)	9100
Gene ID	9100
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human USP10 (Q14694).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14694
Protein Weight	87kDa
Shipping	Ice bag