

Usp10 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14694

Item Number: CM0028047



Introduction

Rabbit polyclonal antibody targeting human USP10, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse. USP10 is a deubiquitinase regulating p53, autophagy, ribosome recycling, and NF-kB signaling. Ideal for cancer and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	USP10 Rabbit pAb
Remarks/Alias	UBPO; USP10
Species	Human
Gene ID (Human)	9100
Gene ID	9100
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 499-798 of human USP10 (NP_005144.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14694
Protein Weight	87kDa
Shipping	Ice bag