

Usp39 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q53Gs9

Item Number: CM0006742



Introduction

Validate USP39 (Q53GS9) with this rabbit monoclonal antibody. Suitable for WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, and rat. Targets a deubiquitinase involved in pre-mRNA splicing, DNA repair, antiviral response, and NF-κB regulation.

[Learn More](#)

Parameter	Value
Product Name	USP39 Rabbit mAb
Remarks/Alias	65K; SAD1; CGI-21; HSPC332; SNRNP65; USP39
Species	Human
GeneID (human)	10713
GeneID	10713
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 466-565 of Human USP39 (Q53GS9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53GS9
Protein Weight	65kDa
Shipping	Ice bag