


Usp7/Hausp Rabbit Pab - Q93009

Item Number: CM0016624



Introduction

USP7/HAUSP rabbit polyclonal antibody for WB, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Targets deubiquitinating enzyme USP7 in p53 regulation and DNA repair. Ideal for protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	USP7/HAUSP Rabbit pAb
Remarks/Alias	TEF1; HAUSP; HAFOUS; DEL16P13.2; C16DELp13.2; USP7/HAUSP
Species	Human
GeneID (human)	7874
GeneID	7874
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-1102 of human USP7/HAUSP (NP_003461.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93009
Protein Weight	128 kDa
Shipping	Ice bag