


Usp9X Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q93008

Item Number: CM0004842



Introduction

USP9X Rabbit Monoclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat USP9X, a deubiquitinase critical for protein turnover and signaling.

[Learn More](#)

Parameter	Value
Product Name	USP9X Rabbit mAb
Remarks/Alias	FAF; FAM; hFAM; DFFRX; MRX99; XLID99; MRXS99F; USP9X
Species	Human
Gene ID	8239
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human USP9X (Q93008).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93008
Protein Weight	290kDa
Shipping	Ice bag