

Uty Rabbit Polyclonal Antibody For Wb, Elisa - O14607

Item Number: CM0026026



Introduction

Rabbit polyclonal antibody against human UTY (histone demethylase UTY), validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat, suitable for epigenetics and Y chromosome research.

[Learn More](#)

Parameter	Value
Product Name	UTY Rabbit pAb
Remarks/Alias	UTY1; KDM6C; KDM6AL; UTY
Species	Human
Gene ID (Human)	7404
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human UTY (NP_009056.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14607
Protein Weight	150kDa
Shipping	Ice bag