

B-Synuclein Rabbit Pab For Wb, Elisa - Q16143

Item Number: CM0017099

Introduction

β -Synuclein Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human β -synuclein, with cross-reactivity to mouse and rat. Ideal for research on Alzheimer's disease and other neurodegenerative disorders.

[Learn More](#)

Parameter	Value
Product Name	β -Synuclein Rabbit pAb
Remarks/Alias	SNCB; β -Synuclein
Species	Human
GeneID (Human)	6620
GeneID	6620
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human β -Synuclein (NP_001001502.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16143
Protein Weight	14kDa
Shipping	Ice bag