

Vamp2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P63027

Item Number: CM0012863



Introduction

High-quality rabbit polyclonal antibody against VAMP2 (synaptobrevin-2), validated for WB, IHC-P, and ELISA. Recognizes human, mouse, rat VAMP2. Ideal for neurotransmitter release and synaptic vesicle fusion studies.

[Learn More](#)

Parameter	Value
Product Name	VAMP2 Rabbit pAb
Remarks/Alias	SYB2; VAMP-2; NEDHAHM; VAMP2
Species	Human
Gene ID (Human)	6844
Gene ID	6844
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VAMP2 (NP_055047.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63027
Protein Weight	13kDa
Shipping	Ice bag