

## Vangl2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9ULK5

Item Number: CM0026080



### Introduction

Rabbit polyclonal antibody targeting VANGL2 (Strabismus 1), validated for WB, IHC-P, ELISA. Detects human and mouse VANGL2. Ideal for planar cell polarity, neural tube development, and morphogenesis research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | VANGL2 Rabbit pAb                                                                                                 |
| Remarks/Alias    | LPP1; LTAP; STB1; STBM; STBM1; VANGL2                                                                             |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 57216                                                                                                             |
| GeneID           | 57216                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human VANGL2 (NP_065068.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse                                                                                                      |
| SWISS            | Q9ULK5                                                                                                            |
| Protein Weight   | 60kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |