

Vars Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P26640

Item Number: CM0025075

Introduction

High-quality polyclonal antibody against human VARS (Valyl-tRNA synthetase) validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. 140 kDa. Ideal for translation and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	VARs Rabbit pAb
Remarks/Alias	G7A; VARS; VARS2; NDMSCA
Species	Human
Gene ID (Human)	7407
Gene ID	7407
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human VARS (NP_006286.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26640
Protein Weight	140kDa
Shipping	Ice bag