

Vcam1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P19320

Item Number: CM0007658



Introduction

High-quality VCAM1 rabbit monoclonal antibody validated for WB, IHC-P, and ELISA. Recognizes human VCAM1 (P19320) with immunogen derived from amino acids 25-698. Ideal for inflammation and immune surveillance research.

[Learn More](#)

Parameter	Value
Product Name	VCAM1 Rabbit mAb
Remarks/Alias	CD106; INCAM-100; VCAM1
Species	Human
Gene ID (Human)	7412
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-698 of human VCAM1 (NP_001069.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P19320
Protein Weight	71KDa/74KDa/81KDa
Shipping	Ice bag