


Vcp Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P55072

Item Number: CM0026010



Introduction

High-quality rabbit polyclonal antibody targeting VCP (p97/TERA) for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat VCP. Ideal for ERAD, autophagy, and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	VCP Rabbit pAb
Remarks/Alias	p97, TERA, CDC48, FTDALS6, VCP
Species	Human
GeneID	7415
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-770 of human VCP (NP_009057.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55072
Protein Weight	89kDa
Shipping	Ice bag