


Vegfc Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P49767

Item Number: CM0017694



Introduction

Anti-VEGFC rabbit pAb validated for WB, IHC-P, ELISA. Recognizes human and mouse VEGFC, a key angiogenic and lymphangiogenic growth factor. Ideal for vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	VEGFC Rabbit pAb
Remarks/Alias	VRP, Flt4-L, LMPH1D, LMPHM4, VEGFC
Species	Human
Gene ID (Human)	7424
Gene ID	7424
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human VEGFC (NP_005420.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49767
Protein Weight	47kDa
Shipping	Ice bag