

Veph1 Rabbit Polyclonal Antibody For Wb And Elisa - Q14D04

Item Number: CM0015871

Introduction

VEPH1 rabbit polyclonal antibody for Western blot and ELISA. Targets human VEPH1 (MELT) protein, cross-reacts with mouse. Suitable for TGF-beta signaling research.

[Learn More](#)

Parameter	Value
Product Name	VEPH1 Rabbit pAb
Remarks/Alias	MELT; VEPH; VEPH1
Species	Human
Gene ID (Human)	79674
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 448-547 of human VEPH1 (NP_078897.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14D04
Protein Weight	95kDa
Shipping	Ice bag