


Bii-Tubulin/ β 2-Tubulin Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13885 / Q9Bva1

Item Number: CM0009206



Introduction

Rabbit monoclonal antibody against β II-Tubulin/ β 2-Tubulin (TUBB2A/TUBB2B) validated for WB, IHC-P, IF/ICC, and ELISA; reacts with human, mouse, rat; ideal for microtubule dynamics research.

[Learn More](#)

Parameter	Value
Product Name	β II-Tubulin/ β 2-Tubulin Rabbit mAb
Remarks/Alias	CDCBM5; TUBB; TUBB2; β II-Tubulin/ β 2-Tubulin
Species	Human
Gene ID (Human)	7280/347733
Gene ID	7280/347733
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 346-445 of human β II-Tubulin/ β 2-Tubulin (Q13885).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13885/Q9BVA1
Protein Weight	50kDa
Shipping	Ice bag