

Villin1 Rabbit Polyclonal Antibody For Wb, Elisa - P09327

Item Number: CM0027871



Introduction

High-quality rabbit polyclonal antibody against human Villin-1, validated for WB and ELISA. Detects endogenous Villin-1 in human, mouse, and rat samples. Ideal for epithelial cell and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Villin1 Rabbit pAb
Remarks/Alias	VIL; D2S1471; Villin1
Species	Human
Gene ID (Human)	7429
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-827 of human Villin1 (NP_009058.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09327
Protein Weight	93kDa
Shipping	Ice bag