

Vimp Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bqe4

Item Number: CM0014427



Introduction

VIMP rabbit polyclonal antibody validated for WB and ELISA, reactive with human, mouse, and rat. Targets selenoprotein S involved in ER-associated degradation. Ideal for endoplasmic reticulum stress research.

[Learn More](#)

Parameter	Value
Product Name	VIMP Rabbit pAb
Remarks/Alias	SELS; VIMP; ADO15; SBB18; SEPS1; AD-015
Species	Human
GeneID (Human)	55829
GeneID	55829
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human VIMP (NP_982298.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BQE4
Protein Weight	21kDa
Shipping	Ice bag