


Vinculin Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P18206

Item Number: CM0017750



Introduction

Vinculin Rabbit Polyclonal Antibody (PAb) validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat vinculin (~124kDa). Suitable for cell adhesion, cytoskeleton, and mechanosensing studies.

[Learn More](#)

Parameter	Value
Product Name	Vinculin Rabbit pAb
Remarks/Alias	MV; MVCL; CMD1W; CMH15; HEL114; Vinculin
Species	Human
Gene ID (Human)	7414
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 751-850 of human Vinculin (NP_054706.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18206
Protein Weight	124kDa
Shipping	Ice bag