


Vip Rabbit Pab For Wb, Ihc-P, Elisa - P01282

Item Number: CM0017706



Introduction

VIP Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat VIP. Ideal for neuropeptide signaling research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	VIP Rabbit pAb
Remarks/Alias	PHM27; VIP
Species	Human
GeneID (human)	7432
GeneID	7432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-170 of human VIP (NP_003372.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01282
Protein Weight	19 kDa
Shipping	Ice bag