


Vip Rabbit Polyclonal Antibody For Wb, Elisa - P01282

Item Number: CM0022524



Introduction

VIP rabbit polyclonal antibody validated for Western blot (WB) and ELISA, suitable for human, mouse, and rat samples. Targets the VIP neuropeptide involved in vasodilation, neuroprotection, and immune modulation.

[Learn More](#)

Parameter	Value
Product Name	VIP Rabbit pAb
Remarks/Alias	PHM27; VIP
Species	Human
GeneID (Human)	7432
GeneID	7432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-169 of human VIP (NP_919416.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01282
Protein Weight	19kDa
Shipping	Ice bag