


Vista Rabbit Pab - Q9H7M9

Item Number: CM0015627



Introduction

High-quality rabbit polyclonal antibody against human VISTA (B7H5, GI24) for WB and ELISA applications. Suitable for T-cell activation research.

[Learn More](#)

Parameter	Value
Product Name	VISTA Rabbit pAb
Remarks/Alias	B7H5; GI24; B7-H5; Dies1; PD-1H; SISP1; VISTA; PP2135; C10orf54; DD1alpha
Species	Human
Gene ID	64115
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 33-194 of human VISTA (NP_071436.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H7M9
Protein Weight	34kDa
Shipping	Ice bag