


Vldlr Rabbit Pab For Wb, Elisa - P98155

Item Number: CM0030111



Introduction

VLDLR Rabbit pAb, specific for human and mouse very low-density lipoprotein receptor, suitable for WB and ELISA. Key reagent for lipid metabolism, Reelin signaling, and neuronal migration research. Sourced from CamelBio's one-stop IVD raw material platform.

[Learn More](#)

Parameter	Value
Product Name	VLDLR Rabbit pAb
Remarks/Alias	CAMRQ1; CARMQ1; CHRMQ1; VLDL-R; VLDLRCH; VLDLR
Species	Human
GeneID (Human)	7436
GeneID	7436
Immunogen	Recombinant protein corresponding to a sequence containing amino acids 833-869 of human VLDLR (NP_003374.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P98155
Protein Weight	96kDa
Shipping	Ice bag