

Vsnl1/Vilip-1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62760

Item Number: CM0028214

Introduction

High-quality VSNL1/VILIP-1 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for neuroscience and Alzheimer's disease biomarker research.

[Learn More](#)

Parameter	Specification
Product Name	VSNL1/VILIP-1 Rabbit pAb
Remarks/Alias	HLP3; VILIP; HPCAL3; HUVISL1; VILIP-1; VSNL1/VILIP-1
Species	Human
GeneID (Human)	7447
GeneID	7447
Immunogen	Recombinant fusion protein containing amino acids 2-191 of human VSNL1/VILIP-1 (NP_003376.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62760
Protein Weight	22kDa
Shipping	Ice bag