

WASF3 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Upy6

Item Number: CM0021186

Introduction

Rabbit polyclonal antibody targeting human WASF3 (WAVE3). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying actin cytoskeleton dynamics and cell morphology. - Q9UPY6

[Learn More](#)

Parameter	Value
Product Name	WASF3 Rabbit pAb
Remarks/Alias	SCAR3; WAVE3; Brush-1; WASF3
Species	Human
Gene ID	10810
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human WASF3 (NP_001278894.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPY6
Protein Weight	55kDa
Shipping	Ice bag