


Wdr77 Mouse Monoclonal Antibody For Wb And Elisa - Q9Bqa1

Item Number: CM0004401



Introduction

Mouse monoclonal antibody against human WDR77 (p44/Mep50), validated for Western blot and ELISA applications. Suitable for investigating methylosome function, spliceosomal Sm protein assembly, and Piwi protein methylation.

[Learn More](#)

Parameter	Value
Product Name	WDR77 Mouse mAb
Remarks/Alias	p44; MEP50; MEP-50; HKMT1069; Nbla10071; p44/Mep50; WDR77
Species	Human
Gene ID (Human)	79084
Gene ID	79084
Immunogen	Recombinant protein of human WDR77
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BQA1
Protein Weight	37kDa
Shipping	Ice bag