

Wisp3 Rabbit Polyclonal Antibody For Wb, Elisa - 095389

Item Number: CM0017206

Introduction

High-quality rabbit polyclonal antibody against human WISP3 (CCN6), validated in WB and ELISA, with cross-reactivity to mouse and rat. Ideal for studying mitochondrial respiration and skeletal growth.

[Learn More](#)

Parameter	Specification
Product Name	WISP3 Rabbit pAb
Remarks/Alias	PPD; LIBC; PPAC; PPRD; WISP3; WISP-3
Species	Human
GeneID (Human)	8838
GeneID	8838
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 42-260 of human WISP3 (NP_937882.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95389
Protein Weight	39kDa
Shipping	Ice bag