


Wnk3 Rabbit Mab - Q9Byp7

Item Number: CM0006493



Introduction

Rabbit monoclonal antibody against human WNK3, validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat WNK3 (~198 kDa). Ideal for kinase signaling and electrolyte homeostasis research.

[Learn More](#)

Parameter	Value
Product Name	WNK3 Rabbit mAb
Remarks/Alias	PRKWNK3; WNK3
Species	Human
Gene ID	65267
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence within amino acids 1700-1800 of human WNK3. (NP_065973.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYP7
Protein Weight	198 kDa
Shipping	Ice bag