

Wtap Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15007

Item Number: CM0020858



Introduction

WTAP Rabbit Polyclonal Antibody (pAb) validated for WB, IF/ICC, ELISA with cross-reactivity to human, mouse, rat. Targets pre-mRNA splicing regulator WTAP, a key m6A methylation factor. Ideal for RNA modification and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	WTAP Rabbit pAb
Remarks/Alias	Mum2; WTAP
Species	Human
GeneID (Human)	9589
GeneID	9589
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 212-355 of human WTAP (NP_004897.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15007
Protein Weight	44kDa
Shipping	Ice bag