


Xbp1 Rabbit Monoclonal Antibody For Wb, Ip, And Elisa - P17861

Item Number: CM0012563



Introduction

XBP1 rabbit monoclonal antibody suitable for Western blot, immunoprecipitation, and ELISA. Detects human and mouse XBP1. Useful for studying ER stress and unfolded protein response.

[Learn More](#)

Parameter	Value
Product Name	XBP1 Rabbit mAb
Remarks/Alias	XBP2; TREB5; XBP-1; TREB-5; XBP1s
Species	Human
GeneID (Human)	7494
GeneID	7494
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17861
Protein Weight	29kDa
Shipping	Ice bag