


Xbp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P17861

Item Number: CM0023633



Introduction

Rabbit polyclonal antibody against human XBP1 (SWISS: P17861). Validated for WB, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Suitable for ER stress, UPR, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	XBP1 Rabbit pAb
Remarks/Alias	XBP2; TREB5; XBP-1; TREB-5; XBP1
Species	Human
GeneID (Human)	7494
GeneID	7494
Immunogen	Recombinant protein A synthetic peptide corresponding to a sequence within amino acids 277-376 of human XBP1 (NP_005071.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17861
Protein Weight	29kDa
Shipping	Ice bag