

Xpo7 Rabbit Pab For Wb And Elisa - Q9Uia9

Item Number: CM0021979

Introduction

XPO7 Rabbit pAb is a polyclonal antibody against human exportin-7 (RANBP16), validated for WB and ELISA. Reacts with human, mouse, and rat. Used for studying nuclear transport mechanisms.

[Learn More](#)

Parameter	Value
Product Name	XPO7 Rabbit pAb
Remarks/Alias	EXP7; RANBP16; XPO7
Species	Human
GeneID (Human)	23039
GeneID	23039
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1000 to the C-terminus of human XPO7 (NP_055839.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UIA9
Protein Weight	124kDa
Shipping	Ice bag