


Xrcc4 Rabbit Polyclonal Antibody For Wb - Q13426

Item Number: CM0020672



Introduction

Rabbit anti-XRCC4 polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and rat XRCC4 (38 kDa). Essential tool for DNA double-strand break repair and NHEJ pathway studies.

[Learn More](#)

Parameter	Value
Product Name	XRCC4 Rabbit pAb
Remarks/Alias	SSMED, hXRCC4, XRCC4
Species	Human
GeneID (Human)	7518
GeneID	7518
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-336 of human XRCC4 (NP_001304941.1?))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13426
Protein Weight	38kDa
Shipping	Ice bag