

Yeats2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ulm3

Item Number: CM0023051



Introduction

YEATS2 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, and rat YEATS2, a chromatin reader that binds H3K27 crotonylation. Essential tool for epigenetics and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	YEATS2 Rabbit pAb
Remarks/Alias	FAME4; YEATS2
Species	Human
Gene ID (Human)	55689
Gene ID	55689
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-331 of human YEATS2 (NP_060493.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9ULM3
Protein Weight	151kDa
Shipping	Ice bag