

## YKL-40/CHI3L1 Rabbit Pab - P36222

Item Number: CM0019026



### Introduction

Anti-YKL-40/CHI3L1 rabbit polyclonal antibody validated for WB, IF/ICC, IHC-P, ELISA applications. Detects human, mouse, rat CHI3L1, a 43 kDa secreted glycoprotein involved in inflammation and tissue remodeling.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | YKL-40/CHI3L1 Rabbit pAb                                                                                                  |
| Remarks/Alias    | GP39; ASRT7; GP-39; YK-40; YKL40; CGP-39; YKL-40; YYL-40; HC-gp39; HCGP-3P; hCGP-39; YKL-40/CHI3L1                        |
| Species          | Human                                                                                                                     |
| GeneID (human)   | 1116                                                                                                                      |
| GeneID           | 1116                                                                                                                      |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-220 of human YKL-40/CHI3L1 (NP_001267.2) |
| Source           | Rabbit                                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                                     |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | P36222                                                                                                                    |
| Protein Weight   | 43 kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                   |