

Ykt6 Rabbit Polyclonal Antibody For Wb And Elisa - O15498

Item Number: CM0020506



Introduction

Rabbit polyclonal antibody targeting human YKT6, validated for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for investigating SNARE-mediated vesicle trafficking and Golgi apparatus organization.

[Learn More](#)

Parameter	Value
Product Name	YKT6 Rabbit pAb
Remarks/Alias	YKT6
Species	Human
GeneID	10652
Immunogen	Synthetic Peptide; A synthetic peptide corresponding to a sequence within amino acids 100-198 of human YKT6 (NP_006546.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15498
Protein Weight	22 kDa
Shipping	Ice bag