


Yod1 Rabbit Polyclonal Antibody For Wb And Elisa - Q5Vvq6

Item Number: CM0014227



Introduction

Validated YOD1 Rabbit pAb for WB and ELISA applications, recognizing human and mouse YOD1 ubiquitin thioesterase OTU1. Involved in ERAD and autophagy pathways. Suitable for IVD research. - Q5VVQ6

[Learn More](#)

Parameter	Value
Product Name	YOD1 Rabbit pAb
Remarks/Alias	DUBA8; OTUD2; PRO0907; YOD1
Species	Human
GeneID (Human)	55432
GeneID	55432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-348 of human YOD1 (NP_061036.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5VVQ6
Protein Weight	38kDa
Shipping	Ice bag