


## Ythdc2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H6S0

Item Number: CM0016242



### Introduction

Rabbit polyclonal antibody against human YTHDC2, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying m6A RNA helicase function in meiosis and spermatogenesis.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | YTHDC2 Rabbit pAb                                                                                                     |
| Remarks/Alias    | CAHL; hYTHDC2; YTHDC2                                                                                                 |
| Species          | Human                                                                                                                 |
| GeneID (human)   | 64848                                                                                                                 |
| GeneID           | 64848                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1131-1430 of human YTHDC2 (NP_073739.3) |
| Source           | Rabbit                                                                                                                |
| Category         | Polyclonal Antibodies                                                                                                 |
| Application      | WB, ELISA                                                                                                             |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | Q9H6S0                                                                                                                |
| Protein Weight   | 160kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                               |