

Ythdf2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y5A9

Item Number: CM0020132



Introduction

Rabbit polyclonal antibody against human YTHDF2, an m6A RNA reader regulating mRNA stability. Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for epitranscriptomics and RNA biology research.

[Learn More](#)

Parameter	Value
Product Name	YTHDF2 Rabbit pAb
Remarks/Alias	DF2, CAHL, HGRG8, NY-REN-2, YTHDF2
Species	Human
Gene ID (Human)	51441
Immunogen	A synthetic peptide corresponding to amino acids 143-243 of human YTHDF2 (NP_001166599.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5A9
Protein Weight	62kDa
Shipping	Ice bag