

Ywhaz/14-3-3 Zeta/Delta Rabbit Pab For Wb, Elisa - P63104

Item Number: CM0018139



Introduction

Rabbit polyclonal antibody targeting human YWHAZ/14-3-3 zeta/delta, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying signaling pathways, phosphorylation, and protein interactions.

[Learn More](#)

Parameter	Value
Product Name	YWHAZ Rabbit pAb
Remarks/Alias	HEL4; YWHAD; KCIP-1; HEL-S-3; POPCHAS; HEL-S-93; 14-3-3-zeta; YWHAZ
Species	Human
GeneID (Human)	7534
GeneID	7534
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63104
Protein Weight	28kDa
Shipping	Ice bag