

Yy1 Rabbit Pab For Wb, Elisa, Chip - P25490

Item Number: CM0018480



Introduction

YY1 rabbit polyclonal antibody for WB, ELISA, ChIP. Reacts with human, mouse, rat YY1 protein (~45 kDa). Ideal for studying transcription regulation, chromatin remodeling, and DNA repair. - P25490

[Learn More](#)

Parameter	Value
Product Name	YY1 Rabbit pAb
Remarks/Alias	DELTA; NF-E1; UCRBP; GADEV5; INO80S; YIN-YANG-1; YY1
Species	Human
GeneID (Human)	7528
GeneID	7528
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 80-160 of human YY1 (NP_003394.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P25490
Protein Weight	45kDa
Shipping	Ice bag