

Zak Rabbit Pab - Zak

Item Number: CM0024915

Introduction

ZAK Rabbit pAb is a polyclonal antibody raised in rabbit against recombinant human ZAK protein. It detects ZAK in human, mouse, and rat samples via WB, IHC-P, IF/ICC, and ELISA, supporting research on MAP3K signaling and stress response.

[Learn More](#)

Parameter	Value
Product Name	ZAK Rabbit pAb
Remarks/Alias	pk; AZK; MLT; MRK; ZAK; CNM6; MLK7; MLTK; SFMMP; mlklak; MLTKbeta; MLTKalpha
Species	Human
Gene ID (Human)	51776
Gene ID	51776
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human ZAK (NP_598407.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
Protein Weight	35kDa/51kDa/91kDa
Shipping	Ice bag